



# ProteinExpress

ProteinExpress Co., Ltd.  
Chiba University Inohana Innovation Plaza  
1-8-15, Inohana, Chuo-ku,  
Chiba-shi, Chiba 260-0856, Japan  
Tel: +81-43-202-5755, Fax: +81-43-202-5756  
E-mail: [antibody@proteinexpress.co.jp](mailto:antibody@proteinexpress.co.jp)  
Web: [www.proteinexpress.co.jp](http://www.proteinexpress.co.jp)

## Anti-Mouse SRPK2 (KIBB3384) Polyclonal Antibody, Rabbit

**Code # :** MKB3384 [FSR07744], (Lot# XXX)

**Gene :** mouse SFRS protein kinase 2 (mSRPK2, mKIBB3384)

**Immunogen :** GX2929 (GST-fusion protein, 170 amino acids)  
NPLDPRNADKIRVKIADLGNACWVHKHFTEDIQTRQYRSIEVLIGAGYSTPADIWSTACMAFELA  
TGDYLFEPHSGEDYSRDEDHIAHIIELLGSIPRHFALS<sup>1</sup>SGKYSREFFNRRGELRHITKLPWSLFDV  
LVEKYGWPHE<sup>2</sup>DAAQFTDFLIPMLEMVPEKRASAGECLRH

**Format :** Rabbit Whole Serum

**Preservative :** 0.05 % NaN<sub>3</sub>

**Quantity :** 100 µL

**Specificity :** Specific to recombinant protein GX2929. This antibody detects mSRPK2 protein. It also recognizes human SRPK2 protein. Other species have not been tested.

**Application :** Western blotting (1 : 1,000), Other applications have not been tested.

**Storage :** Store at -20°C. Avoid freeze-thaw cycles.

**References :** 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

**For research use only, not for diagnostic use**



**Kazusa DNA Research Institute**  
2-6-7, Kazusa-kamatari,  
Kisarazu, Chiba 292-0818 JAPAN  
<http://www.kazusa.or.jp>