



COSMO BIO CO., LTD.  
Inspiration for Life Science

## Anti-KIAA0302 Rat Homologue Gene Product ( $\beta$ SpIII) Polyclonal Antibody, Rabbit

- Code** : PRX-PBR-1003
- Format** : Protein A purified Rabbit IgG
- Presentation** : Lyophilized, reconstitute with 0.2 ml of distilled water
- Immunogen** : GST fusion-KIAA0302 Rat Homologue Gene Product ( $\beta$ SpIII)(292 amino acids: Q2097-K2388):  
QPPTSEPMASQPEGLVDGQRVLDTAWDGTQSKLPPSTQA  
PSINGVCTDTESSQPLLEQQRLEQSNVPEGPGSGTGDDESS  
GPRGERQTLPRGPAPSPMPQSRSSSAHVATLPARGAELS  
AQEQMEGTLCRKQEMEAFNKKAAANRSWQNVYCVLRRGSL  
GFYKDARAASAGVPYHGEVPVSLARAQGSVAFDYRKRKHV  
FKLGLQDGKEYLFQAKDEAEMSSWLRVVNAAIATASSASGE  
PEEPVPSASRGLTRAMTAPPVVSQPEGSIVLRSKDGRERE  
REKRFSFFKKNK
- Buffer** : 0.01 M Phosphate Buffer containing 0.15 M NaCl, pH 7.4
- Preservatives** : 0.05 %  $\text{NaN}_3$
- Quantity** : 0.1 mg / vial
- Applications** : Immunohistochemistry 1:5
- Storage** : 4°C

For research use only, not for diagnostic use

Manufactured by ProteinExpress Co., Ltd.



ProteinExpress



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619