



## Anti PACAP38 (Human) Serum Cat. No. YII-Y040-EX Lot No. 083171013

**Description:** This antiserum was raised in a rabbit by immunization with a carrier free synthetic PACAP38 (hu- man) peptide. The product vial contains 50μL of the titled antiserum obtained by lyophilizing its 0.001 M phos- phate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohistochemistry or any other immu- noreaction with PACAP.

Immunogen: Synthetic PACAP38 (human), carrier free Host: Rabbit

Amino Acid Sequence of PACAP38 (human)<sup>1)</sup>: HSDGIFTDSY SRYRKQMAVK KYLAAVLGKR YKQRVKNK-NH2

Product Form: Lyophilized unpurified serum Size: 50 μL

**Reconstitution:** Reconstitute the product with 0.5mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN3 0.1%).

**Storage:** The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon recon- stitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing- thawing should be avoided.

**Suggested Working Dilution Range:** 1:1,000-2,000 (final dilution ~1:14,000) for radioimmunoassay;

1: 200-1,000 for immunohistochemistry (frozen or paraffin sections)<sup>2)</sup>. Optimal dilution should be determined by each laboratory for each application.

Specificity (based on radioimmunoassay): PACAP38 (human) 100%, PACAP27 (human) 0%, VIP 0%, PHI (rat) 0%, PHI (human) 0%, secretin 0%, helodermin 0%, CCK-8 0%, met-enkephalin 0%, leu-enkaphalin 0%, dynorphin 0%, oxytocin 0%, vasopressin 0%

Positive Control (immunohistochemistry): Rat colon

Species Tested: Rat, ginseng radix<sup>3)</sup>

## REFERENCES:

- 1) A. Miyata, A. Arimura et al., Isolation of a novel 38 residue-hypothalamic polypeptide which stimulates adenylate cyclase in pituitary cells. Biochemistry Biophysics Research Communications 164:567-574, 1989
- 2) S. Kimura, Y. Ohshige et al., Localization of pituitary adenylate cyclase-activating polypeptide (PACAP) in the hypothalamus-pituitary system in rats: light and electron microscopic immunocytochemical studies. Journal of Neuroendocrinology 6: 503-507, 1994
- 3) T. Ohishi, N. Yanaihara et al., Immunoreactivities of neuropeptides in plants (*Ginseng Radix* and *Panacis Japonici Rhizoma*). Proceedings of 18th Gut Hormone Conference, Japan Society of Gut Hormones (Ed) 13:285-291, 1995

## FOR RESEARCH LABORATORY USE ONLY

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

