



## Anti Peptide YY (PYY)(Human) Serum Cat. No. YII-Y072-EX Lot No. 37140301

**Description:** This antiserum was raised in a rabbit by immunization with a carrier free synthetic PYY (human) peptide. The product vial contains 50 μL of the titled antiserum obtained by lyophilizing its 0.001 M phosphate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohistochemistry, or any other immunore- action with PYY.

Immunogen: Synthetic PYY (human), carrier free Host: Rabbit

Amino Acid Sequence of PYY (human)1):

YPIKPEAPGE DASPEELNRY YASLRHYLNL VTRQRY-NH2

Product Form: Lyophilized unpurified serum Size: 50 μL

**Reconstitution:** Reconstitute the product with 0.5mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN3 0.1%).

**Storage:** The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon recon- stitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing- thawing should be avoided.

Suggested Working Dilution Range: 1:3,000-10,000 (final dilution ~1:70,000) for radioimmunoassay;

1: 1,000- 9,000 for immunohistochemistry (frozen or paraffin sections). Optimal dilution should be determined by each laboratory for each application.

**Specificity** (based on radioimmunoassay): PYY (human) 100%, PYY (19-36)(human) <0.001%, PYY (1-19) (human) <0.001%, PYY (rat) 80%, PYY (13-36)(rat) <0.1%, PP (human) 0%, NPY (human) <0.001%

Positive Control (immunohistochemistry): Rat ileum

Species Tested: Rat

## REFERENCES:

1) K. Tatemoto, I. Nakano, et al., Isolation and primary structure of human peptide YY. Biochemistry Biophysics Research Communications 157:713-718, 1988

## FOR RESEARCH LABORATORY USE ONLY

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

