



## Anti GLP-2 (Rat) Serum Cat. No. YII-Y322-EX Lot No. 652150417

**Description:** This antiserum was raised in a rabbit by immunization with a keyhole lympet hemocyanin (KLH) conjugate of synthetic GLP-2 (rat) peptide. The product vial contains  $50 \,\mu\text{L}$  of the titled antiserum obtained by lyo- philizing its 0.001M phosphate buffer (pH 7.0, 0.5mL) solution. It can be used for immunoassay, immunohistoche- mistry or any other immunoreaction with GLP-2 (rat).

Immunogen: Synthetic GLP-2 (rat)-KLH conjugate Host: Rabbit

Amino Acid Sequence of GLP-2 (rat)1):

HADGSFSDEM NTILDNLATR DFINWLIQTK ITD

Product Form: Lyophilized unpurified serum Size: 50 μL

**Reconstitution:** Reconstitute the product with 0.5mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN3 0.1%).

**Storage:** The product will be stable for over one year if it be stored at -20°C to -80°C until opened. Upon recon- stitution, the antiserum solution must be stored at 2°C to 8°C and used within one month. Repeated freezing- thawing should be avoided.

Suggested Working Dilution Range: 1:3,000-50,000 for enzyme immunoassay; 1:500-2,000 for immunohi- stochemistry (frozen or paraffin section). Optimal dilution should be determined by each laboratory for each appli- cation.

**Specificity** (based on enzyme immunoassay): GLP-2 (rat) 100%, GLP-2 (3-33) (rat) 90%, GLP-2 (14-33) (rat) 36%, GLP-2 (hu- man) 100%, GLP-2 (1-18) (rat) < 0.1%, GLP-1 (7-36)-NH2 0%, glucagon 0%, glicentin (rat) 0%

Positive Control (immunohistochemistry): Rat pancreas, rat ileum

Species Tested: Rat

## REFERENCES:

- 1) S. Mojsov, G. Heinrich et al., Preproglucagon gene expression in pancreas and intestine diversifies at the level of post-translational processing. Journal of Biological Chemistry 261:11880-11886, 1986
- 2) I. Kato, J. Li et al., Synthesis of rat glucagon-like peptide (GLP)-2 and its biological and immunochemical studies, Peptide Science, N. Fujii (Ed), The Japanese Peptide Society, p159-162, 2000
- 3) J. Li, N. Yanaihara et al., ELISA for rat and human GLP-2 and the application. Proceeding of 21th Gut Hormone Conference, Japan Society of Gut Hormones (Ed) 63-67, 2001

## FOR RESEARCH LABORATORY USE ONLY

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

