



## Anti Urocortin 2 (Mouse) Serum Cat. No. YII-Y363-EX Lot No. 1344170121

**Description:** This antiserum was raised in a rabbit by immunization with porcine thyroglobulin (pTG) conjugate of synthetic urocortin 2 (mouse) peptide. The product vial contains 50mL of the titled antiserum obtained by lyophilizing its 0.001 M phosphate buffer (pH 7.0, 0.5 mL) solution. It can be used for immunoassay, immunohistochemistry or any other immunoreaction with urocortin 2 (mouse).

Immunogen: Synthetic Urocortin 2 (mouse)-pTG conjugate Host:Rabbit Amino Acid Sequence of Urocortin 2 (mouse)<sup>1</sup>:

VILSLDVPIG LLRILLEQAR YKAARNQAAT NAQILAHV-NH<sub>2</sub>

Product Form: Lyophilized unpurified serum Size: 50mL

**Reconstitution:** Reconstitute the product with 0.5 mL of 0.01M PBS (pH 7.0) to make a 10 fold diluted stock solution. If it is stored in a refrigerator, add moderate antiseptic to the solution (e.g. NaN<sub>3</sub> 0.1%).

**Storage:** The product will be stable for over one year if it be stored at -20 to -80 until opened. Upon reconstitution, the antiserum solution must be stored at 2 to 8 and used within one month. Repeated freezing- thawing should be avoided.

**Suggested Working Dilution Range:**1:8,000 for enzyme immunoassay (EIA); 1:5,000 for immunohistochemistry. Optimal dilution should be determined by each laboratory for each application.

**Specificity**(based on EIA): Urocortin 2 (mouse) 100%, urocortin 2 (7-38) (mouse) 80.52%, urocortin 2 (14-38) (mouse) 24.23%, urocortin 2 (21-38) (mouse) 0.21%, urocortin 2 (26-38) (mouse) 0.41%, urocortin 2 (32-38) (mouse) 0.19%, urocortin 2 (1-13) (mouse) 0.12%, urocortin 1 (mouse, rat) 0%, urocortin 1 (human) 0%, urocortin 2 (rat) < 1%, urocortin 3 (mouse) 0%, ACTH (human) 0%, CRF (human, mouse, rat) 0%, CRF (17-41) (human, mouse, rat) 0%

Positive Control (Immunohistochemistry): Hypothalamus

Species Tested (Immunohistochemistry): Mouse

## Related Peptides:

Urocortin 1 (mouse, rat) DDPPLSIDLT FHLLRTLLEL ARTQSQRERA EQNRIIFDSV-NH<sub>2</sub> Urocortin 3 (mouse) FTLSLDVPTN IMNILFNIDK AKNLRAKAAA NAQLMAQI-NH<sub>2</sub> **Related Antiserum:**Anti Urocortin 1 (mouse, rat) Serum Y362 Anti Urocortin 3 (mouse) Serum Y364

## **REFERENCE**:

1) TM. Reyes, K. Lewis et al., Urocortin II: a member of the corticotropin-releasing factor (CRF) neuropeptide family that is selectively bound by type 2 CRF receptors. *Proc. Natl. Acad. Sci. U. S. A.*, **98**: 2843-2848, 2001

## FOR RESEARCH LABORATORY USE ONLY

DO NOT USE ORGANIC SOLVENTS FOR DISSOLVING ANTISERUM

COSMO BIO CO., LTD. Inspiration for Life Science TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp www.cosmobio.co.jp