

**Monoclonal Antibody to Ia antigen**

Code No.: 7576
Quantity: 200 μ l(200 μ g)
Price: ¥40,000
Clone: ISCR3
Quality: BALB/c IgG2b κ (Purified by protein A chromatography)
Solution: 1mg/ml IgG, 0.1%NaN₃, PBS, pH7.4
Character:

Genetic specificity of interspecies crossreactive ISCR3

Mouse strain	H-2 haplotype	I-region	Reactivity
A.TL	t ₁	k	+
B10.BR	k	k	+
B10.D2	d	d	+
B10.RIII	r	r	+
B10.SM	v	v	-
A.TH	t ₂	s	-
A.SW	s	s	-
B10	b	b	-
B10.M	f	f	-
B10.G	q	q	-

Storage: The solution should be stored below -20°C.

Repeated freezing and thawing should be avoided.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: http://www.cosmobio.co.jp/index_e.asp e-mail: export@cosmobio.co.jp
[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610
FAX : +81-3-5632-9618 FAX : +81-3-5632-9619