



Monoclonal Antibody to **Human Smooth Muscle Myosin Heavy Chain(SM1)**

Abbreviation : MCA MH-01

Code No. : 7599

Price : 200 μ l ······ ¥40,000

Clone : 3F8

Quality : BALB/c IgG1 κ

Solution : Ascites (200 μ l)

Character : Working dilution is about 1:3000

Storage & The solution should be stored below -20°C .

Handling : Repeated freezing and thawing should be avoided.

Note : Research use only.

Reference : Aikawa, M. et al., Circ. Res. 73, 1000-1012 (1993)