



Antibody to Human Atrial Natriuretic Peptide

Abbreviation : RP ANP

Code No. : 7623

Price : 50 μ l ····· ¥50,000

Quality : Rabbit serum

Solution : Serum, Saline (50 μ l)

Character :		Reactivity
to	Human ANP (1-28)	100 %
	Human BNP (1-32)	< 0.01 %
	Human CNP (1-22)	< 0.01 %

Example : Add successively.
 100 μ l of sample in the assay buffer*
 100 μ l of assay buffer
 100 μ l of $\times 1000$ diluted antiserum
 Incubate at 4°C for 24 hours.
 Add 50 μ l of 125 I-HANP solution.
 Incubate at 4°C for 24 hours.
 Separate by the double antibody method.

*Assay buffer

100 mM phosphate buffer pH 7.4
 1 mM EDTA dipotassium
 10 mM glycine
 10 mM ϵ -aminocaproic acid
 1 mM sodium azide
 0.6 % triton X-100

Storage & Handling : The solution should be stored below -20°C .
 Repeated freezing and thawing should be avoided.

Note : Research use only.

