



# Antibody to Arginine Vasopressin

Abbreviation : RP AVP

Code No. : 7626

Price : 50  $\mu$ l..... ¥50,000

Quality : Rabbit serum

Solution : Serum, Saline (50  $\mu$ l)

Character :		Reactivity
	to Arginine Vasopressin	100 %
	Lysine Vasopressin	< 0.3 %
	Oxytocin	< 0.01 %

Example : Add successively.

100  $\mu$ l of sample in the assay buffer\*

100  $\mu$ l of  $\times$ 1000 diluted antiserum

Incubate at 4°C for 24 hours.

Add 100  $\mu$ l of <sup>125</sup>I-AVP solution.

Incubate at 4°C for 24 hours.

Separate by the double antibody method.

\*Assay buffer

100 mM phosphate buffer pH 7.4

10 mM EDTA dipotassium

20 mM glycine

10 mM  $\epsilon$ -aminocaproic acid

1 mM sodium azide

0.1 % triton X-100

0.1 % BSA

Storage & The solution should be stored below -20°C.

Handling : Repeated freezing and thawing should be avoided.

Note : Research use only.

