



Peroxidase Conjugated Monoclonal Antibody to

Human IgG4

Abbreviation: MCA DG-01-HRP  
Code No. : 7634  
Quantity: 50  $\mu$ L (50  $\mu$ g)  
Clone: 5C3  
Quality: BALB/c IgG1 $\kappa$ -Horseradish Peroxidase  
Solution: 1 mg/mL Conjugate, 1% BSA, PBS, pH 7.4

Character:		Reactivity
to	Human IgA	< 0.5%
	Human IgG1	< 0.5%
	Human IgG2	< 0.5%
	Human IgG3	< 0.5%
	Human IgM	< 0.5%

Storage: The solution should be stored below -20°C. Repeated freezing and thawing should be avoided.

Note : Research use only



COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN

[International/Outside Japan]

TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618

e-mail: export@cosmobio.co.jp / URL: www.cosmobio.com

[国内連絡先]

TEL: 03-5632-9610 / FAX: 03-5632-9619

e-mail: support@cosmobio.co.jp / URL: www.cosmobio.co.jp