




Monoclonal Antibody to Human PTH M-region

- Code No.:** 7643
Quantity: 200 μ l (200 μ g)
Price: ¥40,000
Clone: 6D7
Quality: BALB/c mouse IgG₁ κ (Purified by protein A chromatography)
Solution: 1mg/ml IgG, 0.1%NaN₃, PBS, pH7.4
Character: This antibody is specific to mid region (39-68) of human PTH.



- Storage:** The solution should be stored below -20°C.
Repeated freezing and thawing should be avoided.
- Note :** Research use only.