



Monoclonal Antibody to Human PTH C-terminal

Code No.: 7648

Quantity: 200 μ l (200 μ g)


Price: ¥40,000

Clone: 5A7

Quality: BALB/c mouse IgG₁ κ (Purified by protein A chromatography)

Solution: 1mg/ml IgG, 0.1%NaN₃, PBS, pH7.4

Character: This antibody is specific to C-terminal (69-84) of human PTH.



Storage: The solution should be stored below -20°C.

Repeated freezing and thawing should be avoided.

Note: Research use only.