



YAMASA CORPORATION

DIAGNOSTICS DEPARTMENT

Tokyo Office

1-23-8, Nihonbashi-kakigaracho

Chuo-ku, Tokyo, 103-0014, Japan


Tel. +81-3-3668-8558 Fax. +81-3-3668-8407

www.yamasa.com/shindan/index.htm

Monoclonal Antibody to Apolipoprotein B

Code No.: 7649
Quantity: 200 μ L (200 μ g)
Clone: 1D2
Quality: BALB/c IgG1 κ (Purified by protein A chromatography.)
Solution: 1mg/mL IgG, 0.1%NaN₃, PBS, pH7.4
Character :

Specificity of anti-human lipoprotein B 1D2



Storage: The solution should be stored below -20°C.
Repeated freezing and thawing should be avoided.

Note : Research use only.



COSMO BIO Co., LTD.

Inspiration for Life Science