




### Monoclonal Antibody to Apolipoprotein AI

**Code No.:** 7697  
**Quantity:** 200 $\mu$ l(200 $\mu$ g)  
**Price:** ¥40,000  
**Clone:** 2A6  
**Quality:** BALB/c IgG1 $\kappa$  (Purified by protein A chromatography)  
**Solution:** 1mg/ml IgG, 0.1%NaN<sub>3</sub>, PBS, pH7.4  
**Character:**



**Storage:** The solution should be stored below -20°C.

**Repeated freezing and thawing should be avoided.**

**Note ::** Research use only.