



**Peroxidase Conjugated Monoclonal Antibody to  
Mouse Immunoglobulin Kappa Chain**

Code No.: 7962  
 Quantity: 200 $\mu$ l(200 $\mu$ g protein)  
 Price: ¥30,000  
 Clone: 2B6  
 Quality: Horseradish peroxidase conjugated Fisher Rat IgG1.  
 Solution: 1mg/ml protein, 1% Bovine serum albumin, 0.1% ProClin 300 in PBS, pH7.4  
 Storage: The solution should be stored below -20°C.  
 Repeated freezing and thawing should be avoided.  
 Note: Research use only.  
 Character:




Table. Specificities of rat anti-mouse immunoglobulin antibodies

| Clone | Code.<br>Purified/HRP Ab | Test clone <sup>a</sup> / isotype |                       |                       |                        |                       |                        |                       |                        |                      |
|-------|--------------------------|-----------------------------------|-----------------------|-----------------------|------------------------|-----------------------|------------------------|-----------------------|------------------------|----------------------|
|       |                          | 1F4<br>IgG1 $\kappa$              | 2B11<br>IgG1 $\kappa$ | 1D1<br>IgG1 $\lambda$ | 5B11<br>IgG1 $\lambda$ | 3H8<br>IgG2a $\kappa$ | 2G11<br>IgG2a $\kappa$ | 2F2<br>IgG2b $\kappa$ | 2A10<br>IgG2b $\kappa$ | 1D8<br>IgG3 $\kappa$ |
| 14G11 | 7965/7966                | +                                 | +                     | +                     | +                      | -                     | -                      | -                     | -                      | -                    |
| 7D5   | 7967/7968                | -                                 | -                     | -                     | -                      | +                     | +                      | -                     | -                      | -                    |
| 2B6   | 7965/7966                | +                                 | +                     | -                     | -                      | +                     | +                      | +                     | +                      | +                    |
| 7G12  | 7963/7964                | -                                 | -                     | +                     | +                      | -                     | -                      | -                     | -                      | -                    |

<sup>a</sup> Isotype identified mouse monoclonal antibodies were obtained from Yamasa monoclonal antibody libraries. Anti-mouse immunoglobulin antibodies were applied to mouse monoclonal antibodies-coated ELISA plates. Rat anti-mouse immunoglobulin antibodies that reacted the antigens were probed by horseradish peroxidase-labeled anti-rat IgG and TMB, a peroxidase substrate. +, absorbance > 1.2; - absorbance < 0.15.