




**Monoclonal Antibody to  
 Mouse Immunoglobulin G2a**

**Code No.:** 7967  
**Quantity:** 500 µl (500 µg)  
**Price:** ¥40,000  
**Clone:** 7D5  
**Quality:** Fisher rat IgG2a (Purified by protein G chromatography).  
**Solution:** 1mg/ml IgG, 0.1%NaN<sub>3</sub>, PBS, pH7.4  
**Storage:** The solution should be stored below -20°C.  
 Repeated freezing and thawing should be avoided.  
**Note:** Research use only.  
**Character:**



**Table. Specificities of rat anti-mouse immunoglobulin antibodies**

Clone	Code. Purified/HRP Ab	Test clone <sup>a</sup> / isotype								
		1F4 IgG1κ	2B11 IgG1κ	1D1 IgG1λ	5B11 IgG1λ	3H8 IgG2ακ	2G11 IgG2ακ	2F2 IgG2bκ	2A10 IgG2bκ	1D8 IgG3κ
14G11	7965/7966	+	+	+	+	-	-	-	-	-
7D5	7967/7968	-	-	-	-	+	+	-	-	-
2B6	7965/7966	+	+	-	-	+	+	+	+	+
7G12	7963/7964	-	-	+	+	-	-	-	-	-

<sup>a</sup> Isotype identified mouse monoclonal antibodies were obtained from Yamasa monoclonal antibody libraries. Anti-mouse immunoglobulin antibodies were applied to mouse monoclonal antibodies-coated ELISA plates. Rat anti-mouse immunoglobulin antibodies that reacted the antigens were probed by horseradish peroxidase-labeled anti-rat IgG and TMB, a peroxidase substrate. +, absorbance > 1.2; - absorbance < 0.15.

Information

YAMASA CORPORATION  
 Diagnostics Department  
 1-23-8 Nihonbashi Kakigaracho, Chuo-ku, Tokyo 103-0014, Japan  
 Tel: 81-3-3668-8558 Fax: 81-3-3668-8407  
 E-mail: dimac@yamasa.com