




## Peroxidase Conjugated Monoclonal Antibody to Mouse Transferrin

**Code No.:** 7970  
**Quantity:** 200  $\mu$ l (200  $\mu$ g protein)  
**Price:** ¥30,000  
**Clone:** 1H10  
**Quality:** Horseradish peroxidase conjugated fisher rat IgG1 $\kappa$   
**Solution:** 1mg/ml protein, 1% Bovine serum albumin, 0.1% ProClin 300 in PBS, pH7.4  
**Storage:** The solution should be stored below -20°C.  
 Repeated freezing and thawing should be avoided.  
**Note:** Research use only.  
**Character:**

Titration of 7970



Specificity of ELISA using 7619 and 7970 for  
mouse transferrin

