



**Peroxidase Conjugated Monoclonal Antibody to
 Mouse Hemoglobin**

Code No.: 7971
Quantity: 200 µl (200 µg protein)
Price: ¥30,000
Clone: 7A7
Quality: Horseradish peroxidase conjugated fisher rat IgG2aκ
Solution: 1mg/ml protein, 1% Bovine serum albumin, 0.1% ProClin 300 in PBS, pH7.4
Storage: The solution should be stored below -20°C.
 Repeated freezing and thawing should be avoided.
Note: Research use only.
Character:

