



---

## Monoclonal Antibody to Rat Surfactant Protein D

Code No. : 80075  
Product name: MAB TO RAT SURFACTANT PROTEIN D (5F4)  
Quantity: 200  $\mu$ L (200  $\mu$ g)  
Clone: 5F4  
Quality: Mouse IgG2a (Purified by protein A chromatography)  
Solution: 1 mg/mL IgG, PBS, pH 7.4 (200  $\mu$ L)  
Character: This antibody reacts with mouse SP-D.  
It is possible to use for western blotting and ELISA.  
Storage: The solution should be stored below -20°C.  
Repeated freezing and thawing should be avoided.  
Note: Research use only.



COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.com>

e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

Phone : +81-3-5632-9617

FAX : +81-3-5632-9618