



COSMO BIO Co., LTD.
Inspiration for Life Science



YAMASA CORPORATION

DIAGNOSTICS DEPARTMENT

Tokyo Office

1-23-8, Nihonbashi-kakigaracho

Chuo-ku, Tokyo, 103-0014, Japan

Tel. +81-3-3668-8558 Fax. +81-3-3668-8407

www.yamasa.com/shindan/index.htm

Monoclonal Antibody to Rat S100A8

Code No.: 80083

Product name: MAB TO RAT S100A8 (2H6)

Quantity: 200 μ L (200 μ g protein)

Clone: 2H6

Quality: BALB/c IgG1 κ (Purified by Protein A chromatography.)

Solution: 1mg/mL protein in PBS, pH7.4

Character: This antibody reacts to rat S100A8.

It is possible to use for western blotting and immunohistochemistry.

Storage: The solution should be stored below -20°C.

Repeated freezing and thawing should be avoided.

Note: Research use only.

Reference: Koike A, et al., Inflammation.2011 Apr 13.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: http://www.cosmobio.co.jp/index_e.asp E-mail: export@cosmobio.co.jp

[Outside Japan]

[国内連絡先]

Phone : +81-3-5632-9617

Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619