




Recombinant Mouse C-X-C motif chemokine 3 protein (Cxcl3), partial (Active)

Product Code	CSB-AP001091MO
Abbreviation	Recombinant Mouse Cxcl3 protein, partial (Active)
Uniprot No.	Q6W5C0
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered PBS, pH 7.4
Product Type	Chemokine
Immunogen Species	Mus musculus (Mouse)
Biological Activity	Fully biologically active when compared to standard. The biological activity determined by a chemotaxis bioassay using human CXCR2 transfected human 293 cells is in a concentration range of 10-100 ng/ml.
Purity	>97% as determined by SDS-PAGE.
Sequence	AVVASELRCQ CLNTLPRVDF ETIQSLTVTP PGPHCTQTEV IATLKDGQEV CLNPQGPRLQ IIIKKILKSG KSS
Research Area	Immunology
Source	E.coli
Target Names	Cxcl3
Expression Region	28-100aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	7.9 kDa
Protein Length	Partial
PubMed ID	14764687; 16141072; 15489334

Image



**Endotoxin**

Less than 1.0 EU/μg as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.