




# Recombinant Human Thymic stromal lymphopoietin protein (TSLP) (Active)

<b>Product Code</b>	CSB-AP002121HU
<b>Abbreviation</b>	Recombinant Human TSLP protein (Active)
<b>Uniprot No.</b>	Q969D9
<b>Form</b>	Lyophilized powder
<b>Storage Buffer</b>	Lyophilized from a 0.2 µm filtered 20 mM PB, pH 7.4, 150 mM NaCl
<b>Product Type</b>	Colony Stimulating Factor
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Biological Activity</b>	Fully biologically active when compared to standard. The ED50 as determined by a cell proliferation assay using human IL-7Rα and human TSLP R co-transfected murine BaF3 pro-B cells is less than 0.3 ng/ml, corresponding to a specific activity of >3.3x10 <sup>6</sup> IU/mg.
<b>Purity</b>	>98% as determined by SDS-PAGE.
<b>Sequence</b>	M+YDFTNCDFE KIKAAYLSTI SKDLITYMSG TKSTEFNNTV SCSNRPHCLT EIQSLTFNPT AGCASLAKEM FAMKTKAALA IWCPGYSETQ INATQAMKKR RKRVTTNKC LEQVSQLQGL WRRFNRPLLK QQ
<b>Research Area</b>	Immunology
<b>Source</b>	E.coli
<b>Target Names</b>	TSLP
<b>Expression Region</b>	29-159aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag-Free
<b>Mol. Weight</b>	15.1 kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>PubMed ID</b>	11418668; 11480573; 15372022; 15489334; 17242164; 24850429
<b>Image</b>	


**Endotoxin**

Less than 1.0 EU/μg as determined by LAL method.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.