




Recombinant Human Interleukin-36 gamma (IL36G) (Active)

Product Code	CSB-AP004561HU
Abbreviation	Recombinant Human IL36G protein (Active)
Uniprot No.	Q9NZH8
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 μ m Filtered 20 mM Tris-HCl, 100 mM NaCl, 0.1 mM EDTA, pH 8.0
Product Type	Interleukin
Immunogen Species	Homo sapiens (Human)
Biological Activity	The ED50 as determined by its ability to induce IL-8 secretion in A431 human epithelial carcinoma cells is 5-20 ng/ml.
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	SMCKPITGTINDLNQQVWTLQGQNLVAVPRSDSVTPVTVAVITCKYPEALEQG RGDPIYLGIQNPCLYCEKVGEQPTLQLKEQKIMDLYGQPEPVKPFIFYRAK TGRTSTLESVAFPDWFIASSKRDQPIILTSELGKSYNTAFELNIND
Research Area	Immunology
Source	E.coli
Target Names	IL36G
Expression Region	18-169aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	17 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

**Endotoxin**

Less than 1.0 EU/μg as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.