




Recombinant Human Arylacetamide deacetylase (AADAC), partial

Product Code	CSB-BP001014HU
Abbreviation	Recombinant Human AADAC protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P22760
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	PDNVEEPWRMMWINAHLKTIQNLATFVELLGLHHFMDSFKVVGSFDEVPPTS DENVTVTETKFNNILVRVYVPKRKSEALRRGLFYIHGGGWCVGSAALSGYDLL SRWTADRLDAVVVSTNYRLAPKYHFPIQFEDVYNALRWFLRKKVLAKYGVNPE RIGISGDSAGGNLAAAVTQQLDDPDVKIKLKIQSLIYPALQPLDVDLPSYQENS NFLFLSKSLMVRFWSEYFTTDRSLEKAMLSRQHVPVESSHLFKFVNWSSLLPE RFIKGHVYNNPNYGSSELAKKYPGFLDVRAAPLLADDNKLRLPLTYVITCQYD LLRDDGLMYVTRLRNTGVQVTHNHVEDGFHGAFSFLGLKISHRLINQYIEWLK ENL
Research Area	Signal Transduction
Source	Baculovirus
Target Names	AADAC
Protein Names	Recommended name: Arylacetamide deacetylase EC= 3.1.1.3
Expression Region	24-399aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	45.6
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.