




Recombinant Mouse Complement factor I (Cfi)

Product Code	CSB-BP005279MO
Abbreviation	Recombinant Mouse Cfi protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q61129
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	RSPSASDLPQEELVDQKCLLQKYTHRSCNKVFCQPWQRCIEGTCICKLPYQC PRAGTPVCAMNGRSYPTYCHQKSFECLHPEIKFSHNGTCAAEGKFNVSLIYGR TKTEGLVQVKLVDQDERMFICKNSWSMAEANVACVDLGFPLGVRDIQGSFNIS GNLHINDTECLHVHCRGVETSLAECAFTKRRTLSNGLAGVVICYKQDADFPTS LSFQCVNGKHIPQEAKACNGVNDQDQSDCLCKGCRGNASLCKSGVCIPDQ YKCNGEVDCITGEDESRCEEDRQQNIPKGLARSAQGEAEIETEEMLTSGM DNERKRIKSLLPKLSCGVKRNTHTRRKRVIGGKPANVGDPWQVAIKDGQRIT CGGIYIGGCWILTAHCVRPSRAHSYQVWTALLDWLKPNSQLGIQTVKRIVH EKYNGATFQNDIALIEMKMHTGKKECELPNSVPACVPWSPYLFQPNDRCIISG WGRGKDNQKVYSLRWGEVDLIGNCSQFYPPDRYYEKEMQCAGTRDGSIDACK GDSGGPLVCEDINNVTVWGVISWGENCGKPEFPGVYTRVANYFDWISYHVG RSLVSQHNV
Research Area	Epigenetics and Nuclear Signaling
Source	Baculovirus
Target Names	Cfi
Protein Names	Recommended name: Complement factor I EC= 3.4.21.45 Alternative name(s): C3B/C4B inactivator Cleaved into the following 2 chains: 1. Complement factor I heavy chain 2. Complement factor I light chain
Expression Region	19-603aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	69.2
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.