




Recombinant Macaca mulatta (Rhesus macaque) Interleukin-10 (IL10),Baculovirus

Product Code	CSB-BP011580MOW
Abbreviation	Recombinant Rhesus macaque IL10 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P51496
Product Type	Recombinant Proteins
Immunogen Species	Macaca mulatta (Rhesus macaque)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	SPGQGTQSENSCTRFPGNLPMLRLDLRDAFSRVKTFQMKDQLDNILLKESLL EDFKGYLGCGALSEMIQFYLEEVMQAENHDPDIKEHVNSLGENLKTLLRLRLR RCHRFLPCENKSKAVEQVKNAFSKLQEKGVYKAMSEFDIFINYIEAYMTMKIQN
Research Area	Others
Source	Baculovirus
Target Names	IL10
Protein Names	Recommended name: Interleukin-10 Short name= IL-10Alternative name(s): Cytokine synthesis inhibitory factor Short name= CSIF
Expression Region	19-178aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal myc-tagged
Mol. Weight	22.7 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.