



Recombinant Human DNA polymerase theta (POLQ), partial

Product Code	CSB-BP018321HU2
Abbreviation	Recombinant Human POLQ protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O75417
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MNLLRRSGKRRRSESGSDSFSGSGGDSSASPQFLSGSVLSPPPGLGRCLKA AAAGECKPTVPDYERDKLLLANWGLPKAVLEKYHSFGVKKMF EWQAECLLLG QVLEGKNLVYSAPTSAGKTLVAELLILKRVLEMRRKALFILPFVSVAKEKKYYLQ SLFQEVGIKVDGYMGSTSPSRHFSSLDIAVCTIERANGLINRLIEENKMDLLGM VVVDELHMLGDSHRGYLLELLLT KICYITRKSASCQADLASSLSNAVQIVGMSA TLPNLELVASWLNAELYHTDFRPVPLLESVKVGNSIYDSSMKLVREFEPMQV KGDEDHVSLCYETICDNHSVLLFCPSKKWCEKLADIAREFYNLHHQAEGLVK PSECPPVILEQKELLEVMQDLRRLPSGLDSVLQKTVPWGVAFHHAGLTFEERD IIEGAFRQGLIRVLAATSTLSSGVNLPARRVIIRTPIFGGRPLDILTYKQMVGRAG RKGVDTVGESILICKNSEKSKGIALLOGSLKPVRSCLQRREGEEVTGSMIRAILE IIVGGVASTSQDMHTYAACTFLAASMKEGKQGIQRNQESVQLGAIEACVMWLL ENEFIQSTEASDGTGK VYHPHLG SATLSSSLSPADTL DIFADLQ RAMKGFVL ENDLHILYLVTMPFEDWTTIDWYRFFCLWEKLPTSMKRVAELVGVEEGFLARC VKGKVVARTERQHRQMAIHKRFFTSVLVLLD LISEVPLREINQKYGCNRRGQIQSL QQSAAVYAGMITVFSNRLGWHNMELLLSQFQKRLTFGIQRELCDLVRVSL LNA QRARVLYASGFHTVADLARANIVEVEVILKNAVFPK SARKAVDEEEEEAVEERR NMRTIWVTGRKGLTEREAAALIVEEARMILQQDLVEMGVQWNPCALLHSSTCS LTHSESEVKEHTFISQTKSSYKKLTSKNKSNTIFSDSYIKHSPNIVQDLNKSREH TSSFNCNFGNGNQEHQTCSIFR
Research Area	Epigenetics and Nuclear Signaling
Source	Baculovirus
Target Names	POLQ
Expression Region	1-987aa




Notes Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Tag Info N-terminal 10xHis-tagged and C-terminal Flag-tagged

Mol. Weight 112.9 kDa

Protein Length Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.