




Recombinant Arabidopsis thaliana

At1g09870/F21M12_26 (At1g09870)

Product Code	CSB-BP3069DOA
Abbreviation	Recombinant Mouse-ear cress Multiple inositol polyphosphate phosphatase 1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q941B2
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Proteins
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	DQGFDVRHHLSTVTRYSTSKDVTQNLIEGSNVPSECTPIHLNLVARHGT RSPT KKRLRELESLAGRFKELVRDAEARKLP SDKIPGWL GQWKSPWEGKVKG GELI RQGEDELYQLGIRVRERFPSLFEEDYHPDVYTIRATQIPRASASAVA FGMGLFS EKG NLGPGRNRAFAVTSENRASDTKL RFFECCQNYKSYRKAKEPAVDKLKEP VLNKITASVAKRYDLKFTKQDISSLWFLCKQEASLLNVTNQSCELFTPSEVALLE WTD DLEVFLKGYGNSLNYKM GVPLLEDVLHSMEEAIKAREEKLPPGSYEKAR LRF AHAETIVPFSCLLGLFLDGSEFEKIQKEKPLELPPQPPKTRDFRGSTMAPF GGNNILVLYSCPAESSPKYFVQVLHNEHPIAVPGCDGKDFCPLEDFKAKVVTP HLKHAFDNL CNADLNDLKQKPASSKLSILSSWLF GSSSHDTEL
Research Area	Others
Source	Baculovirus
Target Names	At1g09870
Expression Region	20-487aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	55.2 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.