




Recombinant Clostridium acetobutylicum

Acetoacetate decarboxylase (adc)

Product Code	CSB-BP338540DUB
Abbreviation	Recombinant Clostridium acetobutylicum Acetoacetate decarboxylase protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23670
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Clostridium acetobutylicum (strain ATCC 824 / DSM 792 / JCM 1419 / LMG 5710 / VKM B-1787)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MLLVNQSHQGFGNKEHTSKMVSAIVLYVLLAAAAHSAFARTMLKDEVIKQISTPL TSPAFPRGPYKFHNREYFNIVYRTDMDALRKVVPEPLEIDEPLVRFEIMAMHDT SGLGCTESGQAIPVSFNGVKGDYLMHMYLDNEPAIAVGRELSAYPKKLGYPK LFVDSDTLVGTLDYGKLRVATATMGYKHKALDANEAKDQICRPNYMLKIIPNYD GSPRICELINAKITDVTVHEAWTGPTRLQLFDHAMAPLNDLPVKEIVSSSHILADI ILPRAEVIYDYLK
Research Area	Others
Source	Baculovirus
Target Names	adc
Expression Region	1-244aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	33.2 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.