




Recombinant Human Guanine nucleotide-binding protein G(q) subunit alpha (GNAQ)

Product Code	CSB-BP342271HUn8
Abbreviation	Recombinant Human GNAQ protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P50148
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Guanine nucleotide-binding protein G(q) subunit alpha (GNAQ)
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	MTLESIMACCLSEEAKEARRINDEIERQLRRDKRDARRELKLLLLGTGESGKST FIKQMRIIHGSGYSDDEKRGFTKL VYQNI FTAMQAMIRAMDTLKIPYKYEHNKA HAQLVREVDVEKVSAFENPYVDAIKSLWNDPGIQECYDRRREYQLSDSTKYYL NDLDRVADPAYLPTQQDVL RVRVPTTGII EYPFDLQSVIFRMVDVGGQRSERR KWIHCFENVTSIMFLValseyDQVLVESDNENRMEESKALFRTIITYPWFQNSS VILFLNKKDLLEEKIMYSHLVDYFPEYDGPQRDAQAAREFILKMFVDLNPDSDKI IYSHFTCATDTENIRFVFAAVKDTILQLNLKEYNLV
Research Area	Cardiovascular
Source	Baculovirus
Target Names	GNAQ
Expression Region	1-359aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal Flag-tagged and C-terminal 8xHis-tagged
Mol. Weight	44.2 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.