




Recombinant Vesicular stomatitis Indiana virus Matrix protein (M)

Product Code	CSB-BP356390VBGd8
Abbreviation	Recombinant Vesicular stomatitis Indiana virus Matrix protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04876
Product Type	Recombinant Proteins
Immunogen Species	Vesicular stomatitis Indiana virus (strain Glasgow) (VSIV)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MSSLKKILGLKGKGKSKKLGIAPPPYEEDTSMEYAPSAPIDKSYFGVDEMDTH DPNQLRYEKSFFTVKMTVRSNRPFRITYSDVAAAVSHWDHMYIGMAGKRPFYK ILAF LGSSNLKATPAVLADQGQPEYHAHCEGRAYLPHRMGKTPPMLNVPEHF RRPFNIGLYKGTIELTMTIYDDESLEAAPMIWDHFNSSKFSDFREKALMFGLIVE EEASGAWVLDSVRHSKWASLASSF
Research Area	Signal Transduction
Source	Baculovirus
Target Names	M
Expression Region	71.4 kDa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal MBP-tagged and C-terminal 6xHis-tagged
Mol. Weight	71.4 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.