




# Recombinant Marmota monax Tumor necrosis factor (TNF), partial

<b>Product Code</b>	CSB-BP523593MQG
<b>Relevance</b>	Cytokine that binds to TNFRSF1A/TNFR1 and TNFRSF1B/TNFR. It is mainly secreted by macrophages and can induce cell death of certain tumor cell lines. It is potent pyrogen causing fever by direct action or by stimulation of interleukin-1 secretion and is implicated in the induction of cachexia, Under certain conditions it can stimulate cell proliferation and induce cell differentiation. The TNF intracellular domain (ICD) form induces IL12 production in dendritic cells.
<b>Abbreviation</b>	Recombinant Marmota monax TNF protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O35734
<b>Alias</b>	Cachectin TNF-alpha Tumor necrosis factor ligand superfamily member 2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Marmota monax (Woodchuck)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	GPQREEFLNNLPLSPQAQMLTLRSSSQNMNDKPVAVHVVAKNEDKEQLVWLS RRANALLANGMELIDNQLVVPANGLYLVYSQVLFGKGQCPSYVLLTHTVSRFA VSYQDKVNLLSAIKSPCPKESLEGAEFKPWYEPIYLGGVFELQKGDRLSAEVN LPSYLDFAESGQVYFGVIAL
<b>Research Area</b>	Immunology
<b>Source</b>	Baculovirus
<b>Target Names</b>	TNF
<b>Expression Region</b>	57-233aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged
<b>Mol. Weight</b>	22.2kDa
<b>Protein Length</b>	Extracellular Domain
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.