




Recombinant Human Protein mono-ADP-ribosyltransferase TIPARP (TIPARP)

Product Code	CSB-BP801800HU
Abbreviation	Recombinant Human TIPARP protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z3E1
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MEMETTEPEPDCVVQPPSPDDFSCQMRLSEKITPLKTCFKKKDQKRLGTGT LRSLRPILNTLLESGSLDGVFRSRNQSTDENSLHEPMMKKAMEINSSCPAEN NMSVLIPDRTNVGDQIPEAHPSTEAPERVVPIQDHSFPSETLSGTVADSTPAHF QTDLLHPVSSDVPTSPDCLDKVIDYVPGIFQENSFTIQYILDTSDKLSTELFQDK SEEASLDLVFELVNQLQYHTHQENGIEICMDFLQGTCIYGRDCLKHHTVLPYH WQIKRTTTQKWQSVFNDSQEHLERFYCNPENDRMRMKYGGQEFWADLNAM NVYETTEFDQLRRLSTPPSSNVNSIYHTVWKFFCRDHFHWREYPESVIRLIEEA NSRGLKEVRFMWNNHYILHNSFFRREIKRRPLFRSCFILLPYLQTLGGVPTQ APPPLEATSSSQIICPDGVTSANFY PETWVYMHPSQDFIQVPVSAEDKSYRIIY NLFHKTVPFEKYRILQILRVQNQFLWEKYKRKKEYMNRKMFGDRDRIINERHLFH GTSQDVVDGICKHNFDPRVCGKHATMFGQGSYFAKKASYSHNFSKKSSKGV HFMFLAKVLTGRYTMGSHGMRRPPPVNPGSVTSDLYDSCVDNFFEPQIFVIFN DDQSYFYFVIQYEEVSNTVSI
Research Area	Epigenetics and Nuclear Signaling
Source	Baculovirus
Target Names	TIPARP
Expression Region	1-657aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	79.0 kDa
Protein Length	Full Length



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.