




Recombinant Human Pleckstrin homology domain-containing family F member 1 (PLEKHF1)

Product Code	CSB-BP842759HU
Abbreviation	Recombinant Human PLEKHF1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96S99
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MVDHLANTEINSQRIAAVESCFGASGQPLALPGRVLLGEGVLTKECRKKAKPRI FFLFNDILVYGSIVLNKRKYRSQHIIPLEEVTLLELLPETLQAKNRWMIKTAKKSFV VSAASATERQEWISHIEECVRRQLRATGRPPSTEHAAAPWIPDKATDICTMRCTQ TRFSALTRRHHCRCGFFVCAECRQRFLLPRLSPKPVRVCSLCYRELAAQQ RQEEAEEQGAGSPGQPAHLARPICGASSGDDDDSDDEDKEGSRDGDWPSSV EFYASGVAWSAFHS
Research Area	Others
Source	Baculovirus
Target Names	PLEKHF1
Protein Names	Recommended name: Pleckstrin homology domain-containing family F member 1 Short name= PH domain-containing family F member 1 Alternative name(s): Lysosome-associated apoptosis-inducing protein containing PH and FYVE domains Short name=
Expression Region	1-279aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	35.2 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.