




Recombinant Human Serine/arginine repetitive matrix protein 2 (SRRM2), partial

Product Code	CSB-BP892369HU
Abbreviation	Recombinant Human SRRM2 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UQ35
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	RTARRGSRSSPEPKTKSRTPPRRRSSRSSPELTRKARLSRRSRSSASSSPETR SRTPPRHRRSPSVSSPEPAEKSRSSRRRRSASSPRTKTTSSRRGRSPSPKPRG LQSRSRSRREKTRTTRRRDRSGSSQSTSRRRQSRSRSRVTRRRRGGS YHSRSPARQESSRTSSRRRRGRSRTPTSRKRSRSTSPAPWKRSRSRASP ATHRRSRRTPLISRRRSRSTSPVSRRRSRSTSVTRRRRSRSPVSRRR SRSTPPVTRRRSRSTPTTRRRSRSTPPVTRRRSRSTPPVTRRRSRST SPITRRRSRSTSPVTRRRSRSTSPVTRRRSRSTSPVTRRRSRSTPPAIR RRSRSTPLLPRKRSRSPLAIRRRSRSTPRTARGKRSLTRSPPAIRRSA SGSSSDRSR
Research Area	Others
Source	Baculovirus
Target Names	SRRM2
Protein Names	Recommended name: Serine/arginine repetitive matrix protein 2 Alternative name(s): 300 kDa nuclear matrix antigen Serine/arginine-rich splicing factor-related nuclear matrix protein of 300 kDa Short name= SR-related nuclear matrix prote
Expression Region	1666-2089aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	53.8 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.