



Recombinant Human Acetolactate synthase-like protein (ILVBL)


Product Code	CSB-CF011687HU
Abbreviation	Recombinant Human ILVBL protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A1L0T0
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose
Product Type	Transmembrane Proteins
Immunogen Species	Homo sapiens (Human)
Sequence	METPAAAAPAGSLFPSFLLLACGTLVAALLGAAHRLGLFYQLLHKVDKASVRH GGENVAA VLRAHGVRFIFTLVGGHISPLLVACEKLGIRVVDTRHEVTAVFAADAMARLSGT VGVAAV TAGPGLTNTVTAVKNAQMAQSPILLGGAASTLLQNRGALQAVDQLSLFRPLC KFCVSVR RVRDIVPTLRAAMAAQSGTPGPVFVELPVDVLYPYFMVQKEMVPAKPPKGL VGRVVSWE LENYLANLFAGAWEPQPEGPLPLDIPQASPQQVQRCVEILSRAKRPLMVLGSQ ALLTPTS ADKLRAAVETLGVPCFLGGMARGLLGRNHPLHIRENRSALKKADVIVLAGTV CDFRLSY GRVLSHSSKIIIVNRNREEMLLNSDIFWKPQEAVQGDVGSFVLKLVEGLQGQT WAPDWVE ELREADRQKEQTFREKAAMPVAQHLNPVQVLQLVEETLPDNSILVVDGGDFV GTAAHLVQ PRGPLRWLDPGAFTLGVGAGFALGAKLCRPDAEVWCLFGDGAFGYSLIEFD TFVRHKIP VMALVGNDAGWTQISREQVPSLGSNVACGLAYTDYHKAAMGLGARGLLLSRE NEDQVVKV LHDAQQQCRDGHVPVVNIIIGRTDFRDGSIIV
Source	in vitro E.coli expression system
Target Names	ILVBL
Protein Names	Recommended name: Acetolactate synthase-like protein EC= 2.2.1.- Alternative name(s): IlvB-like protein
Expression Region	1-632
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged



Protein Length

Full length protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.