




Recombinant Pongo abelii NADPH oxidase 4 (NOX4)

Product Code	CSB-CF015961PYX
Abbreviation	Recombinant Pongo abelii NOX4 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R5C5
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MAVSWRSWLANEGVKHLCLFIWLSMNVLLFWKTFLLYNQGPEYHYLHQMLGL GLCLSRASASVLNLCNCSLILLPMCRTLLAYLRGSQKVPSRRTRRLLDKSRTFHI TCGVTICIFSGVHVAHLVNALNFSVNYSEDFVELNAARYRDEDPRKLLFTTVP GLTGVCMMVVVLFMITASTYAIRVSNYDIFWYTHNLFFVFYMLLTLHVSGGLLKY QTNLDTHTPPGCISLNRRTSSQNISLPEYFSEHFHEPFPEGFSKPEEFTQNTFVKI CMEEPRFQANFPQTWLWISGPLCLYCAERLYRYIRSNKPVTIISVISHPSDVMEI RMVKENFKARPGQYITLHCPVSVALENHPFTLTMCPTETKATFGVHLKIVGDW TERFRDLLPPSSQDSEILPFIQSRNYPKLYIDGPFPGSPFEESLNYEVS LCVAG GIGVTPFASILNTLLDDWKPYKLRRLYFIWVCRDIQSFRWFADLLCMLHNKFW QENRPDYVNIQLYLSQTDGIQKIIGEKYHALNSRLFIGRPRWKLLFDEIAKYNRG KTVGVFCCGPNSLSKTLHKLSNQINSYGTRFEYNKESFS
Research Area	Cardiovascular
Source	in vitro E.coli expression system
Target Names	NOX4
Expression Region	1-578aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	69.7 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.