




Recombinant Human Transient receptor potential cation channel subfamily V member 1 (TRPV1), partial

Product Code	CSB-CF822774HU1
Abbreviation	Recombinant Human TRPV1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NER1
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Human Transient receptor potential cation channel subfamily V member 1 (TRPV1), partial
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	QLSTAKSRTLFGKGDSEEAFPVDCPHEEGELDSCPTITVSPVITIQRPGDGPT GARLLSQDSVAASTEKTLRLYDRRSIFEAVAQNNCQDLESLLLFLQKSKKHLTD NEFKDPETGKTCLLKAMLNLDGQNTTIPLLEIARQTDSLKELVNASYTDSYY KGQTALHIAIERNMALVTLLVENGADVQAAAHGDDFFKKTGRPGFYFGELPL SLAACTNQLGIVKFLQNSWQTADISARDSVGNTVLHALVEVADNTADNTKFV TSMYNEILMLGAKLHPTLKLEELTNKKGMTPLALAAGTGKIGVLAYILQREIQEP ECRHLSRKFTWAYGPVHSSLYDLSCIDTCEKNSVLEVIAYSSSETPNRHDML LVEPLNRLQLQDKWDRFVKRIFYFNFLVYCLYMIIFTMAAYYRPVDGLPPFKMEK TGDYFRVTGEILSVLGGVYFFFRGIQYFLQRRPSMKTFLVDSYSEMLFFLQSLF MLATVVLYFSLHKEYVASMVFSALGWTNMLYYTRGFQQMGIYAVMIEKMILR DLCRFMFVYIVFLFGFSTAVVTIEDGKNDSLPESTSHRWGRGPACRPPDSSY NSLYSTCLELFKFTIGMGDLEFTENYDFKAVFIILLAYVILTYILLNMLIALMGET VNKIAQESKNIWKLQRAITILDTEKSFLKCMRKAFRSGKLLQVGYPDGKDDYR WCFRVDEVNWTWNTNVGIINEDPGNCEGVKRTLFSLSRSSRVSGRHWKNF ALVPLLREASARDRQSAQPEEVYLRQFSGSLKPEDAIEVFKSPAASGEK
Research Area	Neuroscience
Source	in vitro E.coli expression system
Target Names	TRPV1
Expression Region	39-839aa



Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	92.5 kDa
Protein Length	Partial

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-----------------------	---

Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
-------------------	--