




Recombinant Human RNA exonuclease 4 (REXO4)

Product Code	CSB-CF863931HUa6
Abbreviation	Recombinant Human REXO4 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9GZR2
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MGKAKVPASKRAPSSPVAKPGPVKTLTRKKNKKKKRFWKSKAREVSKKPASG PGAVVRPPKAPEDFSQNWKALQEWLLKQKSQAPEKPLVISQMGSKKKPKIIQQ NKKETSPQVKGEEMPAGKDQEASRGVPSGSKMDRRAPVPRTKASGTEHNK KGTKERTNGDIVPERGDIEHKKRKAKEAAPAPPTTEEDIWFDDVDPADIEAAIGP EAAKIARKQLGQSEGSVSLSLVKEQAFGGLTRALALDCEMVGVGPKGEESMA ARVSIVNQYGKCVYDKYVKPTEPVTDYRTAVSGIRPENLKQGEELEVQKEVA EMLKGRILVGHALHNDLKVLFLDHPKKKIRDTQKYKPFKSQVKSGRPSLRLLSE KILGLQVQQA EHC SIQDAQAAMRLYVMVKKEWESMARDRRPLL TAPDHCSDD A
Research Area	Epigenetics and Nuclear Signaling
Source	in vitro E.coli expression system
Target Names	REXO4
Protein Names	Exonuclease XPMC2 Prevents mitotic catastrophe 2 protein homolog Short name: hPMC2
Expression Region	1-422aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-B2M-tagged
Mol. Weight	60.7kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.