



Recombinant Rat Transient receptor potential cation channel subfamily V member 2 (Trpv2)


Product Code	CSB-CF893416RA
Abbreviation	Recombinant Rat Trpv2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9WUD2
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MTSASSPPAFRLETSDGDEEGNAEVNKGKQEPPPMESPFQREDRNSSPQIKV NLNFIKRPPKNTSAPSQQEPDRFDRDLFSVVS RGVPEELTGLLEYLRWNSKY LTDSAYTEGSTGKTCLMKAVLNLQDGVNACIMPLLQIDKDSGNPKLLVNAQCT DEFYQGHSA LHIAIEKRSLQCVKLLVENGADVHLRACGRFFQKHQGTCTFYFGE LPLSLAACTKQWDVVITYLLENPHQPASLEATDSLGN TVLHALVMIADNSPENS ALVIHMYDGLLQMGARLCPTVQLEEISNHQGLTPLKLAAKEGKIEIFRHILQREF SGPYQPLSRKFTEWCYGPVRVSLYDLSSVDSWEKNSVLEIIAFHCKSPNRHR MVLLEPLNKLLQEKWDRLVSRFFNFACYLVYMFIFTVVAYHQPSLDQPAIPSS KATFGESMLLLGHILILLGGIYLLLGQLWYFWRRRLF IWISFMDSYFEILFLLQAL LTVLSQVLRFMETEWYLP LLVLSVLGWLNLLYYTRGFQHTGIYSVMIQKVILR DLLRFLLVYLVFLFGFAVALVSLSREARSPKAPEDNNSTVTEQPTVGQEEEPAP YRSILDASLELFKFTIGMGELAFQEQLRFRGVVLLLLLAYVLLTYVLLLNMLIALM SETVNHVADNSWSIWKLQKAISVLEMENGYWWCRRKKHREGRLKVGTRGD GTPDERWCFRVEEVNWA AWEKTLPTLSEDPSGPGITGNKKNPTSKPGKNSA SEEDHLPLQVLQSP
Research Area	Neuroscience
Source	in vitro E.coli expression system
Target Names	Trpv2
Expression Region	1-761aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged



Mol. Weight 88.2 kDa

Protein Length Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.