




Recombinant Drosophila melanogaster Gustatory receptor for sugar taste 43a (Gr43a)

Product Code	CSB-CF895567DLU
Abbreviation	Recombinant Drosophila melanogaster Gr43a protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9V4K2
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Drosophila melanogaster Gustatory receptor for sugar taste 43a (Gr43a)
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sensitivity	Not Test
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	<p> MEISQPSIGIFYISKVLALAPYATVRNSKGRVEIGRSWLFTVYSATLTVVMVFLT YRGLLFDANSEIPVRMKSATSKVVTALDVSVVVMAIVSGVYCGLFSLNDTLELN DRLNKIDNTLNAYNNFRRDRWRALGMAAVSLLAISILVGLDVGTWMRIAQDMNI AQSDTELVHWHYIPFYSLYFILTGLQVNIANTAYGLGRRFGRLNRMLSSSFLAE NNATSAIKPQKVSTVKNVSVNRPAMPSPALHASLTKLNGETLPSEAAAKNKGLLL KSLADSHESLGKCVHLLSNSFGIAVLFILVSCLLHLVATAYFLFLELLSKRDNGY LWVQMLWICFHFLRLLMVVEPCHLAARESRTIQIVCEIERKVHEPILAEAVKKF WQQLLVVDADFSACGLCRVNRTILTSFASAIATYLVILIQFQRTNG </p>
Research Area	Others
Source	in vitro E.coli expression system
Target Names	Gr43a
Expression Region	1-427aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	48.5 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.