






# Recombinant Mouse Allograft inflammatory factor 1 (Aif1)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                       |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP001490MO                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Relevance</b>         | Actin-binding protein that enhances mbrane ruffling and RAC activation. Enhances the actin-bundling activity of LCP1. Binds calcium. Plays a role in RAC signaling and in phagocytosis. May play a role in macrophage activation and function. Promotes the proliferation of vascular smooth muscle cells and of T-lymphocytes. Enhances lymphocyte migration. Plays a role in vascular inflammation. |
| <b>Abbreviation</b>      | Recombinant Mouse Aif1 protein                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                       |
| <b>Uniprot No.</b>       | O70200                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Alias</b>             | Ionized calcium-binding adapter molecule 1                                                                                                                                                                                                                                                                                                                                                            |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Immunogen Species</b> | Mus musculus (Mouse)                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Purity</b>            | Greater than 90% as determined by SDS-PAGE.                                                                                                                                                                                                                                                                                                                                                           |
| <b>Sequence</b>          | SQSRDLQGGKAFGLLKAQQEERLEGINKQFLDDPKYSNDEDLPSKLEAFKVKY<br>MEFDLNGNGDIDIMSLKRMLEKLGVPKTHLELKRILIREVSSGSEETFSYSDFLR<br>MMLGKRSAILRMILMYEEKNKEHKRPTGPPAKKAISELP                                                                                                                                                                                                                                           |
| <b>Research Area</b>     | Others                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Target Names</b>      | Aif1                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Expression Region</b> | 2-147aa                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                   |
| <b>Tag Info</b>          | N-terminal 6xHis-tagged                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Mol. Weight</b>       | 20.8kDa                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | Full Length of Mature Protein                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Image</b>             |                                                                                                                                                                                                                                                                                                                                                                                                       |



We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.