




Recombinant Rat Allograft inflammatory factor 1 (Aif1)

Product Code	CSB-EP001490RA
Relevance	Actin-binding protein that enhances mbrane ruffling and RAC activation. Enhances the actin-bundling activity of LCP1. Binds calcium. Plays a role in RAC signaling and in phagocytosis. May play an role in macrophage activation and function. Promotes the proliferation of vascular smooth muscle cells and of T-lymphocytes. Enhances lymphocyte migration. Plays a role in vascular inflammation .
Abbreviation	Recombinant Rat Aif1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P55009
Alias	Ionized calcium-binding adapter molecule 1;MRF-1;Microglia response factor
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	SQSKDLQGGKAFLGKLAQQEERLDGINKHFLDDPKYSSDEDLQSKLEAFKTKY MEFDLNGNGDIDIMSLKRMLEKLGVPKTHLELKKLIREVSSGSEETFSYSDFLR MMLGKRSAILRMILMYEEKNKEHQKPTGPPAKKAISELP
Research Area	Neuroscience
Source	E.coli
Target Names	Aif1
Protein Names	Recommended name: Allograft inflammatory factor 1 Short name= AIF-1 Alternative name(s): Ionized calcium-binding adapter molecule 1 MRF-1 Microglia response factor
Expression Region	2-147aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	32.7kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.