




Recombinant Human Aldo-keto reductase family 1 member B10 (AKR1B10) (K125L)

Product Code	CSB-EP001540HU(M)
Abbreviation	Recombinant Human AKR1B10 protein (K125L)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
Uniprot No.	O60218
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Human Aldo-keto reductase family 1 member B10 (AKR1B10) (K125L)
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MATFVELSTKAKMPIVGLGTWKSPLGKVKEAVKVAIDAGYRHIDCAYVYQNEH EVGEAIQEIKIQEKAVKREDLFIVSKLWPTFFERPLVRKA FEKTLKDLKLSYLDVY LIHWPQGFGKSGDDLFPLDDKGNAIGGKATFLDAWEAMEELVDEGLVKALGVS NFSHFQIEKLLNKPGLKYKPVNTQVECHPYLTQEKL IQYCHSKGITVTAYSPLG SPDRPWAKPEDPSLLEDPKIKEIAAKHKKTAAQVLIRFHIQRNVIVIPKSVTPARI VENIQVFDFKLSDEEMATILSFNRNWRACNVLQSSHLEDYPFNAEY
Research Area	Signal Transduction
Source	E.coli
Target Names	AKR1B10
Expression Region	1-316aa(K125L)
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	62.7 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.