




Recombinant Cynops pyrrhogaster Annexin A5

Product Code	CSB-EP001846EQC
Relevance	Calcium/phospholipid-binding protein which promotes mbrane fusion and is involved in exocytosis.
Abbreviation	Recombinant Cynops pyrrhogaster Annexin A5 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P70075
Product Type	Recombinant Protein
Immunogen Species	Cynops pyrrhogaster (Japanese fire-bellied newt)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MACLKGA KGT VQDAPDFNDKEDAETLRHAMKGLGTDEDTILKLLISRSNKQRQ QIALTYKTLFGRDLTDDLKSEL SGKFETLLVALMVP AHLYDACELRNAIKGLGTL ENVIIIEIMASRTAAEVKNIKETYKKEFDSLEKDIVGDTSGNFERLLVSLVQANR DPVGKVDEGQVENDAKALFDAGENK WGTDEETFISILSTRGVGHLRKVFDQY MTISGYQIEESIQSETGGHFEKLLLAVVKSIRSIIQGYLAEVLYNSMKGAGTDDQ TLIRVLVSRSEIDLFNIRQTFRKH YGKSLHAMIQS DTS GDYRNALLLLCGEIDD
Research Area	Others
Source	E.coli
Expression Region	1-323aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	40.0 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.