




Recombinant Human Apolipoprotein C-IV (APOC4)

Product Code	CSB-EP001934HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P55056
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	CQPEAQEGTLSPPPKLKMSRWSLVRGRMKELLETVVNRTRDQWQWFWSPS TFRGFMQTYDDHLRDLGPLTKAWFLESKDSLLKKTHSLCPRLVCGDKDQG
Research Area	Others
Source	E.coli
Target Names	APOC4
Protein Names	Recommended name: Apolipoprotein C-IV Short name= Apo-CIV Short name= ApoC-IVAlternative name(s): Apolipoprotein C4
Expression Region	27-127aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-B2M-tagged
Mol. Weight	25.8
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.