




# Recombinant Human CMRF35-like molecule 7 (CD300LB), partial

<b>Product Code</b>	CSB-EP004919HUb0
<b>Abbreviation</b>	Recombinant Human CD300LB protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A8K4G0
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	IQGPESVRAPEQGSLTVQCHYKQGWETYIKWWCRGVWRWDTCILIIETRGSEQ GEKSDRVSIKDNQKDRTFTVTMEGLRRDDADVYWCGIERRGPDLTQVKVIV DPEGAASTTASSPTNSNMAVFIGSHKRNH
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	CD300LB
<b>Expression Region</b>	18-151aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged
<b>Mol. Weight</b>	18.8 kDa
<b>Protein Length</b>	Partial

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



## **Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.