




Recombinant Human T-cell antigen CD7 (CD7), partial

Product Code	CSB-EP004953HU
Relevance	Not yet known.
Abbreviation	Recombinant Human CD7 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09564
Alias	GP40T-cell leukemia antigen;T-cell surface antigen Leu-9TP41; CD7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	AQEVQQSPHCTTVPVGASVNITCSTSGGLRGIYLRQLGPQPQDIYYEDGVVP TTDRRFRGRIDFSGSQDNLTITMHLQLSDTGTYTCAITEVNVYGSGLTLVLVT EEQSQGWHRCS DAPPRASALPAPPTGSALPDPQTASALPDPPAASALP
Research Area	Immunology
Source	E.coli
Target Names	CD7
Expression Region	26-180aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	32.4kDa
Protein Length	Extracellular Domain

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.