




Recombinant Human Cyclin-dependent kinase 2-associated protein 1 (CDK2AP1)

Product Code	CSB-EP005062HU
Relevance	specific inhibitor of the cell-cycle kinase CDK2.
Abbreviation	Recombinant Human CDK2AP1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O14519
Alias	Deleted in oral cancer 1 ;DOC-1;Putative oral cancer suppressor
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MSYKPNLAAHMPAAALNAAGSVHSPSTSMATSSQYRQLLSYDYGPPSLGYTQG TGNSQVPQSKYAELLAIIIEELGKEIRPTYAGSKSAMERLKRGIHARGLVRECLA ETERNARS
Research Area	Cell Cycle
Source	E.coli
Target Names	CDK2AP1
Expression Region	1-115aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	39.4kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.